

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

00280 - Cutdown Pad Order code: 4104.00280



Parameters Surgery - filter Quantitative unit

Basic Surgical Skills ks

A multi-layered pad containing filled veins used for practising the cutdown procedure.

Skills

- Palpation to locate vein position
- Incising tissue layers
- Inserting a cannula
- Vein injection both distally and proximally

• Tying off veins

Features

- Layered pad simulating the epidermis and subdermal layers
- Three straight veins (each 150mm long) filled with mock blood

Note

It is recommended that the veins are tied off each time so that the next student may repeat the procedure without losing the contents from the veins.