

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082

500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

# Varioshake VS 15 T - orbital-rocking shaker

Order code: 1304.L003059



## Information about product price on demand

Device type

Manufacturer

Maximum speed

Note

Dimension of moving base

Max. amplitude

**Rotational Speed** 

Shaking motion

Maximum load

Quantitative unit

Shaker

LAUDA

50 rpm

Digitální zobrazení

450 × 450 mm

3°

2 - 50 ot./min

3D lurched

15 kg

ks

### **Specifications and features**

- Silent operating orbital-rocking shaker, long-lasting
- User-friendly microprocessor technique
- Gentle and uniform motion
- Excellently applicable for continuous operation
- Digital shaking frequency readout and setting via LC-display, in steps of 0.5 rpm, gentle start-up
- Remaining time indication through LC display
- RS 232 serial interface for trouble-free data transfer
- Clearly laid out control panel for easy operation
- An a.c. motor with overload protection drives the unit
- Compact, low-wear counterbalanced drive mechanism, ensuring dependability for continuous operation at three-dimensional motion
- Constant speed independent of load
- Housing made of high impact strength polystyrene
- Comprehensive range of accessories

### **Applications**

In AIDS research, for diagnostic tests with test plates, or for growing cell cultures and micro-organisms this unit has become indispensable. The unit is well-suited for use in laboratories, incubation and moderating rooms in ambient temperatures between +10 and +50 °C.

### **Optional Accessories:**

3965, 3966, 3967, 3968, 3983, 3984, 3985, 3986, 3987, 3988, 3989