

## **AC8 - Adult Intubation Trainer**

Order code: **4101.AC8**



Information about product price on demand

### Parameters

Quantitative unit

ks

- This trainer has been developed to learn how to perform intubation under lifelike conditions
- If during intubation, too much pressure is applied through a laryngoscope to the upper incisors, an audio signal is heard
- The correct location of the tube, inserted after intubation, may be confirmed by audible respiratory sounds through a stethoscope (not supplied) lateral to the left or right wall of the thorax
- Includes separate model of the larynx for identification of anatomical structures
- Supplied in a carrying case