

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

## RCT basic - Magnetic stirrer with heating Order code: 1308.3810000



Information about product price on demand

Parameters	
Device type	Magnetic stirrers
Speed [rpm]	0/50 - 1500
Max. volume [I]	20
Plate diameter [mm]	135
Heating capacity [W]	600
Heating	YES
Maximum speed	1500 rpm
Number of positions	1
Quantitative unit	ks

The new improved version of the laboratory bestseller with revolutionary technology and new functions:

- new: Integrated temperature control
- new: Incl. PT 1000 temperature sensor (PT 1000.60)
- new: Exact temperature and speed setting via digital display, even when switched off
- new: Set safety temperature limit displayed digitally
- new: Hot Top indicator >> hot surface warning to prevent burns!
- new: Digital error code display
- With adjustable safety circuit of heating plate temperature (50 360 °C)
- Safety magnetic stirrer with heating, suitable for unsupervised operation

- Bushing according to DIN 12878 for connecting a contact thermometer, e.g. ETS-D5, enables precise temperature control
- High level of safety thanks to improved heat control technology
- Enclosed assembly (IP 42) guarantees long service life
- Highly polished aluminium heating plate for optimum heat transfer
- High magnetic adhesion