



Cena bez DPH

1.529,00 Eur

Price with VAT

1.850,09 Eur

Parameters

Quantitative unit

ks

- Rotation and flame sterilization happen automatically. No need for cooling down of inoculating loops in agar.
- When needed - at the push of a button - all 5 inoculating loops can be flame sterilized in one full rotational cycle. The time for flame sterilization can continuously be set from 2 to 10 seconds. Suited for right- and left-handed use.
- Can be used with all standard gas safety-burners (adapters are optionally available).

- Flame sterilization time can be minimized to save gas and to increase lifetime of inoculating loops
- No accidental usage of hot loops. Maintenance-free switches of the flame sterilization carousel are covered inside the device, no malfunction due to dirt or stray light.
- We recommend schuett manufactured inoculating loops made of special stainless steel wire or platinum-iridium wire (90/10). Length 60 mm, diameter of wire 0,5 mm.