

HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

## schuett colonyQuant - Automated colony counter

Order code: 3002.3082002



Information about product price on demand

## **Parameters**

Quantitative unit ks

- The economic solution for the evaluation of Petri dishes. Thanks to its high efficiency we can work with samples continuously. Only
  the press of one button can show us the sample and at the same time, his analysis.
   Switching between evaluation agar and alimentary dishes can be done simply by clicking a mouse.
- Software for the detection of bacteria according to color, size and shape.
- Chamber for placing the sample with backlight and a color CCD-firewire camera.
- During few seconds will be done the evaluation together with displaying the sample and table with values. Sample pictures can be storing in digital form.
- Agar dishes, alimentary dishes, screw cap, analysis of mixed culture transition range.
- 8 colours together on Petri dish.
- Partial evaluation, 90% adjustable in accordance to GLP
- schuett colonyQuant includes: coloured CCD-firewire camera with autozoom and autofocus, connecting firewire cable, line connector, chamber for backlight samples, lighting from bellow or from above, lucid counting software compatible with Windows XP