



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec
Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: <http://www.helago-cz.cz>

LF01008 - Intradermal Injection Simulator

Order code: **4107.LF01008**



Cena bez DPH

391,00 Eur

Price with VAT

473,11 Eur

Parameters

Cannula, injection, puncturing - filter

Horní končetiny

Quantitative unit

ks

Now your students can practice intradermal injections safely and conveniently with the Life/form Intradermal Injection Simulator. The unit consists of a casting molded from a real forearm from the wrist to just below the elbow. Life/form vinyl skin provides realistic feel and appearance to ensure a realistic training experience. The arm features eight sites for practicing intracutaneous injections. If fluid is properly injected, a characteristic skin welt will form. The welt is removed by withdrawing the fluid after practice. Each site is reusable by dozens of students. Includes arm section with eight injection sites, sealant, syringe, teaching guide, and storage box. The included sealant (LF01009) will extend the life of skin; however, skin is not replaceable and the simulator will eventually need replacement. Size: 13" x 11" x 5" (33 x 28 x 13 cm).