

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

FDL115 - Safety drying chamber for limited quantities of solvent

Order code: **2103.90100292**



Cena bez DPH

10.133,00 Eur 12.260,93 Eur

Price with VAT

Parameters

Type of convection
Brand, Manufacturer

Volume [I]

External dimensions $(w \times h \times d)$ [mm]

Quantitative unit

safety Binder 115

830 × 805 × 685

ks

- Temperature range: room temperature plus 5 °C to 300 °C
- APT.line™ preheating chamber technology
- Controller with time-segment and real-time programming
- 2 chrome-plated racks
- Fresh-air monitoring with audible and visual alarm and automatic deactivation of heating
- Replaceable filter cartridge, Class M6
- Class 2 independent adjustable temperature safety device (DIN 12880) with visual alarm
- Computer interface: RS 422

A BINDER heating oven with a silicone- and dust-free inner chamber and a symmetrical air flow provides the perfect environment for standard-compliant drying of paints and coating materials containing solvents. The FDL Series guarantees short drying time with high temperature accuracy and defined ventilation. Every BINDER heating oven thus represents the basis for the highest-quality results and reproducible tests in the fields of chemistry and surface technology. The replaceable fresh air cartridges are easy to change and ensure constant fresh air quality. The safety concept meets the DIN EN 1539-2015 standard.

Drying chambers for paints and other materials containing solvents: BINDER safety drying chambers provide the ideal environment for drying paints and other coating materials containing solvents. There is absolute temperature accuracy in the interior, which is free from silicone and dust, along with homogeneous temperature distribution, regardless of the amount or size of the sample. BINDER drying chambers guarantee quick and even drying, within a wide temperature range of between 5 °C and 350 °C. The drying chambers operate with absolute temperature accuracy, even with high ventilation, thereby providing the best possible basis for precise results and reproducible tests. High levels of operational safety are ensured by the electronic monitoring system and fresh air cartridge, which can be replaced.

Safety drying chambers with maximum operational safety: BINDER drying chambers are reliable and operate with high temperature accuracy. The units feature 60-mm-thick insulation, keeping heat dissipation to a minimum, and a specimen temperature display with a temperature sensor. Thanks to the fresh air monitoring system with automatic shut-down, the units comply with all of the requirements of EN 1539 and they also correspond to protection class IP 33. The tight housing door with 2-point door closure also helps take care of safety. In practice, the drying chambers are extremely versatile owing to their functionality and options for individual adjustment. The safety chambers can be used for plenty of tasks within chemical and surface engineering applications. Alongside being flexible in their uses and accurate, BINDER drying chambers are also really easy to use in practical terms. BINDER drying chambers provide the ideal conditions for norm-compliant drying of substances containing solvents, with maximum safety guaranteed.