

M118-1 Musclemate II

Order code: **4109.M1181**



Information about product price on demand

Parameters

Cannula, injection, puncturing - filter

Horní končetiny

Quantitative unit

ks

Intramuscular and subcutaneous injection simulator.

Practice:

- Preparation – easy to wear simulator, Velcro closure
- Palpation – check the injection site, acromion is included
- Pinch – the muscle is raised by pinching, it helps to find the right injection site
- Intramuscular injection – simulator can analyse not only the injection site, but also its depth; the surface of the simulator is very similar to real skin and the needle marks do not remain
- Subcutaneous injection – it is possible to practise the subcutaneous injection

Features:

- Built-in lamps and a buzzer mark the mistakes made during the injection
 - Correct injection – the green light turns on, the injection site is well chosen and the depth is right
 - Incorrect injection site – if the injection site is badly chosen, the red light blinks and the buzzer rings
 - Incorrect injection depth – if the injection site is well chosen, but the injection is too deep, the red and the green lights blink and the buzzer rings
- The simulator is light and it fits well; thanks to this fact the practice is very realistic
- Use only water!
- Weight: 0,5 kg
- Case dimensions (w × d × h): 32 × 12 × 11 cm
- Accessories: drain bag, strap