

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 F-mail: info@helago-cz.cz

E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

AD-4212A-PT - Pipette Accuracy Tester

Order code: 0102.AD4212APT



Information about product price on demand

Parameters

Capacity (maximum) 110 g

Min. weighing value 0,001 g

Quantitative unit ks

- Automatic mass-to-volume conversion and accuracy judgment using the special "WinCT-Pipette" software
- Judgment by your preferred criteria manufacturer specifications, ISO8655, or your own
- Includes an evaporation trap to minimize errors due to evaporation
- Three models covering a wide volumetric range
- Portable carrying case for testing easily on-site (not recommended when a 0.001 mg resolution is required)
- Printing function to have paper-based test records
- Standard liquid thermometer for precise mass-to-volume conversion
- Standard calibration weight and tweezers to maintain the balance accuracy
- Pipetting training capability
- RS-232C (with a Serial / USB converter)