



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional
Court in Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: <http://www.helago-cz.cz>

Prestan Professional Adult Jaw Thrust - CPR- AED Training Manikin with Jaw Thrust and CPR Monitor

Order code: **4115.PPJT100MMS**



Cena bez DPH

263,00 Eur

Price with VAT

318,23 Eur

Parameters

Construction

Adult torso

Quantitative unit

ks

The PRESTAN Professional Jaw Thrust Head gives instructors the option to teach with the jaw thrust maneuver or head-tilt/chin-lift allowing the student to grasp the manikin's lower jaw and lift with both hands displacing the jaw of the manikin forward and upward. The jaw thrust maneuver is the safest approach to opening the airway of a victim who has a possible head, neck or spinal injury because it can be

accomplished by minimizing neck movement.

Jaw Thrust Heads may be purchased separately and easily fit onto the torso of the original PRESTAN Professional Adult Manikin.

Fully compliant with current industry guidelines, including the most recent AHA Integrated Feedback Directive.









Available in medium or dark skin tones, with or without CPR Feedback and in Single or 4-Pack configurations.


Additional Details:

All PRESTAN Professional Manikins are available with a revolutionary CPR Rate Monitor that allows for instant feedback to both instructors and students. The visual feedback from the CPR Rate Monitor gives students a real life feel of delivering 100-120 compressions per minute. Lights continually signal rate progression, confirming compressions are within the 100- 120 cpm range recommended by the Industry Guidelines. The yellow light will blink if the 120 cpm upper limit is exceeded, warning students to slow down. Instructors can monitor several students quickly and easily.

In addition to the visual CPR Rate Monitor, students will also hear a clicker sound as the chest is pushed to the appropriate depth of 2.0 -2.4 inches allowing them to experience the true force needed to deliver real life chest compressions.

The CPR Rate Monitor may be purchased separately to upgrade any non-monitored PRESTAN Professional Manikin.

Sleep Mode/No Power	
1st Full Compression (Wake mode/LED Test Mode)	
CPM* < 60	
CPM* 60 – 79	
CPM* 80 – 99	
CPM* 100 - 120	
CPM* > 120	
Compressions Stopped	

*CPM = Compressions per minute
 = Flashing

Includes:

- 1 Manikin
- 1 CPR monitor (with two AA batteries)
- 10 Lung Bags
- Instruction Sheet
- 1 American Red Cross-Branded nylon carrying case for 1 manikin