


HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

M99 - Arterial Puncture Wrist Order code: 4103.M99



Parameters

Cannula, injection, puncturing - filter Quantitative unit

Horní končetiny

ks

Radial artery puncture is a common approach for blood collection and artery catheterization; however, it is a challenging skill for operators. This innovative simulator is designed to provide training in artery puncture with true-to-life feeling.

Features:

- Arterial pulsation is palpable.
- Realistic resistance of tissue and artery wall felt with the injection needle.
- Natural flashback of artificial blood into the needle can be observed.
- Puncture skin, pad and artery tube are replaceable.
- No puncture trace of the injection needle remains on the injection site.
- One-touch, leak free connections.
- Easy clean up.

Set includes:

- 1 wrist and hand model
- 1 circulation pump
- 2 sets of puncture skin (inner and outer skin)
- 4 artery tubes
- 1 silicon tape
- 1 jar of coloring powder (red)
- 1 air bulb
- 1 spoon
- 2 plastic jars
- 2 connection tubes
- 1 vinyl sheet
- 1 carrying case