

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

LF01400 - Life/form® Newborn Nursing Skills and ALS Simulator Order code: 4107.LF01400

<image>



Cena bez DPH Price with VAT

Parameters Construction Quantitative unit 1.265,00 Eur 1.530,65 Eur

ks

Newborn manikin

Use for meeting objectives of the Neonatal Resuscitation Program's course curriculum

• Newborn size: 19" L (crown to heel), 7-1/2-lb. weight, 14-1/2" head circumference

 Includes simulator, replaceable airway (2), bilateral lung, collapsed lung, functional (patent) umbilicus, omphalocele, myelomeningocele, hand skins (2), foot skins (2), vein system for hand, vein system for foot, 500 ml fluid supply bag (2), 3 cc syringe with needle, 25 gauge infusion butterfly needle, 25 ml syringe, artificial blood powder, lubricant, baby powder, and hard carry case

Airway, Breathing, and Ventilation

- Anatomical landmarks
- Auscultation of lung sounds during ventilation
- Bilateral lung expansion with visually dramatic chest rise
- ET tube insertion
- Manually control rate and depth of respiration
- Positive pressure ventilation
- Realistic chest compressions
- Sellick maneuver
- Unilateral chest movement option (collapsed lung and mainstem intubation)

Birth Anomalies

- Neural tube defect (Myelomeningocele) easily attaches to the back
- Omphalocele

Chest Tube

• Practice chest tube insertion and care of the chest tube area

Cleaning and Diaper Changing

- Practice cleaning, bathing, and diaper changing
- Waterproof and completely immersible

Gastrointestinal (GI)

• Accepts a OG/NG tube for practicing tube feeding and suctioning

Monitoring

• Various monitors, sensors, and electrodes may be attached and connected

Observation and Measurement

• Observe and measure each part of the body (head, umbilicus, skin, etc.)

Scenarios

• Build your own high-risk delivery, emergency, transport, and NICU scenarios

Skin

- Delicate appearance and texture of real skin
- Soft, durable lifelike skin

Suction

• Practice suctioning nose, oral cavity, and stomach

Umbilicus

- Can be attached and detached easily
- Made of soft, realistic material
- Practice blood withdrawal or infuse fluids
- Patent umbilicus with venous and arterial access

Venous Access

- Accepts IV catheter or butterfly infusion needle
- IV hand and foot with realistic flashback from pressurized system
- Soft, lifelike skins with small lumen tubing for realistic practice