

DL510H - Ultrasound bath DL 510 H

Order code: **2007.7183**



Information about product price on demand

Parameters

Quantitative unit

ks

Description

- Cleaning of sensitive materials
- Mixing of samples
- Homogenizing of samples for analysis
- Fast degassing of solvents for HPLC

Features

- Power control from 20 to 100 %
- Heating from 20 to 80 °C, thermostatically adjustable, LED display
- Optical LED signal in case of excess temperature
- Automatic short-time extra ultrasound during heating up period to avoid delay in boiling - switch-off possible
- Timer for countdown or continuous operation
- DEGAS - quick degassing
- SweepTec® - permanent sound field oscillation ensures even and gentle cleaning
- Automatic safety shut-down 12 hours after last keypress
- Marking of filling level for safe dosage
- Drain valve for simple emptying, one-piece welded drain bow
- Protection class IP 33 at highest possible level through liquid-tight mains cable and housing front without openings

For more information about this product please download product sheet under "Files for download".