



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: <http://www.helago-cz.cz>

70400 - Clean Bleed™ Mat

Order code: **4104.70400**



Cena bez DPH

746,00 Eur

Price with VAT

902,66 Eur

Parameters

Quantitative unit

ks

The brand NEW Clean Bleed™ Mat has been designed to support all simulation scenarios, with any brand of trainer where external blood or fluid loss is a key indicator and measure of the emergency being simulated.

This easy to setup and mess free product has the ability to create the impression of hemorrhages, such as postpartum, post-operative or gunshot trauma, through a controlled flow into a sealed mat. The flow into the mat is adjustable and can be easily managed, creating a realistic scenario where the trainee or clinician is able to identify fluid output, and respond accordingly.

Versatile and efficient, the Clean Bleed™ Mat can be used under or besides simulated patients, full body manikins or task trainers, in both an outdoor and indoor setting. After simulation, a two-way pump system ensures that the fluid can be pumped back into the reservoir provided, and the internal mat can be removed and rinsed for a hassle free simulation experience.

OVERVIEW

- Controllable blood flow (flow rates of between 50ml and 600ml per minute)
- The flow of blood can be measured either visually, or by weighing
- Easy set-up process

- Provides mess-free simulation training
- Suitable for use with simulated patients, task trainers and full body manikins
- Quick release bag allows trainees to replace the bag quickly, as they would in a real-life scenario
- Absorbent pad behaves as a chux pad would
- Housed in a carry bag

REALISM

- Unique sealed blood flow, heightening the realism for simulated scenarios
- High flows: vaginal trauma (including PPH or meconium fluid); abdominal, chest and neck trauma; loss of limb and post-operative bleeds
- Low flows: head or peripheral trauma; post-operative bleeds; meconium fluid

VERSATILITY

- Can be used in conjunction with SP, task trainers or full-body manikins
- Suitable for both indoor and outdoor settings
- Can be used to simulate a range of fluids
- Portable and durable for use both indoors and outdoors
- Bag allows for the use of standard sized chux pads (60cm x 90cm)

SAFETY

- Latex free

SKILLS GAINED

- Visual inspection or body assessment
- Identification and assessment of the speed and color of the bleed or flow
- Assessment of the speed and relevance of trainee interventions
- Quantifying blood loss (QBL)
- Ability to initiate rapid transfusion protocol if hemorrhage is suspected
- Use of PPE to protect against bodily fluids
- Assessment of the ongoing patient condition and impact of prior interventions on the patient

CONTAINS

- 2x Clean Bleed™ Bag
- 10x Clean Bleed™ Absorbency Pad (70402)
- 1x Clean Bleed™ Connectors Set (70403)
- Clean Bleed™ Fluid Reservoir (70404)
- Clean Bleed™ Carry Case (70405)
- Clean Bleed™ Pump (70408)
- Concentrated Venous Blood - Refill (00021)