

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

Phone: 495 220 229, 495 220 394

LF04301 - TERi - Geriatric Patient Skills Trainer

Order code: 4107.LF04301



Information about product price on demand

Parameters

Nursing simulators - filter

Adult

Quantitative unit

ks

A comprehensive male/female elderly patient care trainer for medical simulation.

#### General

- Excellent range of movement
- Enhanced aesthetics
- Light approx. 23 kg
- Pupils (normal, narrowed, dilated)
- Male and female (wig for female)
- Soft male (circumcised and uncircumcised features) and female genitalia

# Skin care

- Normal and Cancerous Skin Moles
- Stage 1 Sacral Ulcer
- Reddened Skin Folds
- Bandaging and Dressing

#### **Movement Assistance and Fall Prevention**

# **Bathing**

- Clothing Changes
- Bed Bathing

## **Daily Living Assistance**

- Eye Irrigation & Optic Drops
- Denture Placement/Removal
- Ear Canal Rinsing
- Hearing Aid Placement/Removal
- Hair Care

## **Nursing Care**

- Ostomy Care
- Urinary Catheter Care
- Stomach
- IV Arm
- Oxygen Administration
- Left BP Arm (optional)

## **Advanced Features**

- Injection Sites Option (IM Arm, Thigh, Buttocks on the Right Side, Subcutaneous injection)
- Gastric procedures (lavage, gastric tube, probe location)
- Nasogastric Tube Procedures (Placement, Lavage, Gavage, Feeding and Suctioning)
- Tracheostomy Care, Airway and Trachea Suction (Lavage and Suction)
- Urinary Catheterization
- Enema (Female)
- Cervical Swab
- Examination of the Pelvis (Female) 3 Different Cervixes
- Prostate Examination (Male) 5 Different Prostates
- Right IV Arm (optional)
- Inspection and Diagnosis of Feet, Care for Wounds and Swelling

- CPR with Performance Quality Metrics
- Carotid Pulse Manual
- Articulated Jaw for Simulation of Lower Jaw Thrust
- Intubation
- Ventilation