

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

BT-CSPA - Arterial Puncture Training Arm Model

ks

Order code: 4119.BTCSPA



4.400,00 Eur Cena bez DPH Price with VAT 5.324,00 Eur

Parameters

Cannula, injection, puncturing - filter Upper limbs

Ouantitative unit

General Features:

- Arm model and skin pad is similar to real human arm
- · Palpate life like pulse at radial artery
- Changeable pulse intensity and rate
- Check the backward flowing of blood while arterial puncture training
- Realistic resistance at the arterial wall during arterial puncture training
- Hard to find needle marks after punctures
- Skin&Vein All-in-one pad; Easy installation and replacement of skin pad
- The color of palm is changing while Allen test
- Easy clean, Easy maintenance; Self-cleaning system

Arterial Puncture Training Arm Model

Arm model and skin pad is similar to real human arm

- Paplate life-like pulse at radial artery
- Check the backward flowing of blood while arterial puncture training
- Realistic resistance at the arterial wall during arterial puncture training

Easy installation and replacement of skin pad

• All-in-one Skin and Artery

Blood Supply

Supply the arterial blood flow similar to human

• Changeable pulse intensity and rate

Practice anywhere, anytime

• User can select AC power or Battery mode

Easy clean, Easy maintenance

Clinical Skill

- The color of palm is changed while the Allen test
- Practice arterial blood sampling
- Practice arterial cannulation (The arterial BP is different from the actual value)

The Specifications and Components

Adult arm manikin

- 1 Arm manikin
- 2 Skin pads

Blood Supply Device

Etc. Components

- 1 Blood supply bottle
- 1 Simulated blood
- 1 Power adapter
- 1 Storage case