

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

Gynaecological Simulator Order code: 4113.BNDGN001



Information about product price on demand

Parameters

Obstetrics, gynecology - filter

Birthing simulators Quantitative unit ks

The MODEL-med® Gynaecological Simulator is a highly realistic simulation of the female reproductive organs. Made from a very life-like material, it allows the medical practitioner to practice everyday procedures such as bi-manual palpation, pap-smears, and the use of the

The MODEL-med® Gynaecological Simulator is so easy to use!

To change over a uterus, simply put a small amount of lubricant on the cervix and slide the uterus into the cavity at the back of the model until it seats down into the vagina. It is that easy. To swap over another uterus, simply slide the current one out and slide in the one you

The MODEL-med® Gynaecological Simulator is made of a very soft, highly elastic, and very durable material. It is easy to clean and following simple maintenance instructions ensures the models will last a long time. They are also easy to repair if they are accidentally damaged.