

## **Jaw Crusher BB 250 jaws steel for heavy-metal free grinding/ wearing plates steel for heavy-metal free grinding, 400V**

Order code: **1306.210010004**



Information about product price on demand

### Parameters

Jaws	steel for heavy-metal free grinding
Wearing plates	steel for heavy-metal free grinding
Accessories	NO
Grind	Jaw crushers
Quantitative unit	ks

The Jaw Crusher BB 250 is used for the rapid, effective crushing and pre-crushing of medium-hard, hard, brittle and tough materials. Its variety of materials offered, including heavy-metal free steel, its efficiency and safety makes this jaw crusher ideal for sample preparation in laboratories and industrial plants.

**Material feed size:** < 120 x 90 mm

**Final fineness:** < 2 mm high crushing ratio

- front door facilitates cleaning

- infinite/continuous gap width setting
- overload protection
- wide range of materials for contamination free grinding
- no-rebound feed hopper, removable
- sample collector with safety switch
- collecting receptacle with outlet for continuous operation
- connector for dust extraction
- special version with automated sorting of undersize (3 fractions) and oversize (1 fraction) (2 - 6 mm available)
- special version Jaw Crusher and Sample Divider Combination Unit ABP 250 available (max. 5.6 l, divided in 8 equal parts or 1 part sample)