

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

## **ROCKER 3D basic, 3D shaker**

Order code: 1305.0400000



Cena bez DPH 1.148,00 Eur
Price with VAT 1.389,08 Eur

Parameters

Device type Mini shaker Maximum speed 30 rpm Max. amplitude 30 mm **Rotational Speed** fixed 30 Shaking motion 3D Timer NO Maximum load 2 kg Quantitative unit ks

Shaker with three-dimensional tumbling motion and a fixed speed for smooth mixing tasks. Used for homogeneous mixing in flasks, culture flasks, Petri dishes and tubes. Ideal for e.g. the cultivation of cell cultures, DANN extractions, distributing of cells and staining and destaining of cells. Even with load and voltage fluctuations, shaking speed remains constant. A fixed angle ensures reproducible test results.

- Different attachments enable use for a variety of applications
- Usable at 4 to 50 °C in incubators
- Anti-slip rubber mat included
- Easy operation with touch keypad
- Suitable for continuous operation
- Shaking diameter: 8° (angle)
- Speed range: 30 rpm (fixed)

## **Scope of delivery**

- ROCKER 3D basic
- RD 1 universal mat