


HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

Trayster basic, overhead shaker

Order code: 1305.04005000



Cena bez DPH 1.079,00 Eur
Price with VAT 1.305,59 Eur

Parameters

Device type Test tube shaker

ManufacturerIKAMaximum speed20 rpmRotational Speedfixed 20Shaking motionFlippingTimerNOMaximum load2 kgQuantitative unitks

Overhead shaker with vertical rotation for smooth but effective mixing with a fixed speed. Ideal for biological samples, particularly blood samples. Suitable for separating and mixing processes of powdered and liquid samples in Eppendorf tubes and Greiner tubes up to 50 ml. The easy removable attachments allow simple and quick cleaning of the unit of spilled samples.

- Different attachments allow a variety of applications
- Usable at 4 to 50 °C in incubators

- Suitable for continuous operation
- Easy operation by foil keypad
- Speed range: 20 rpm (fixed)