



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in  
Hradec Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: [info@helago-cz.cz](mailto:info@helago-cz.cz)  
Web: <http://www.helago-cz.cz>

---

**R66600 - Epidural-Spinal Injection  
Simulator GENESIS**  
Order code: **4003.R66600**



Information about product price on demand

Parameters

Cannula, injection, puncturing - filter

Quantitative unit

Trunk

ks

The GENESIS Epidural-Spinal Injection Simulator is the ideal training simulator for common neuraxial blocks. This ultrasound compatible model contains anatomical structures found in the lumbar region and realistic dural puncture similar to a human patient. With the realistic loss of resistance and fluid filled intrathecal space users are able to locate the correct target sites. GENESIS is the beginning of lumbar block proficiency.

**Features:**

- Spinal space (Fluid filled tubing simulates dural click and CSF confirmation)
- Epidural space (Simulate LOR technique)
- Vertebral bodies (Locate interspace with visible and palpable spinous processes)
- Ligamentum flavum (Tough and gritty texture provides realistic tactile feel)
- Iliac crest (Standard landmark for L4 spinous process)

**Training Techniques:**

- Epidural LOR (Identify Epidural Loss of Resistance)
- Lumbar Puncture (Confirm Spinal Placement with Fluid Indication)
- Catheter Placement (Insert Epidural & Spinal Catheters)
- CSE Technique (Ideal for Combined Spinal Epidural Technique)
- Ultrasound Compatible (Increase Ultrasound Competency)

Additional puncture pad for geriatric and obese training available.

Supplied with two puncture pads and sic spinal tubings