

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

AirSim Difficult Airway Order code: 4110.DA91100



Cena bez DPH 4.860,00 Eur

Price with VAT 5.880,60 Eur

Parameters

Quantitative unit ks

Adult intubation head trainer provides a solution for varying degrees of difficulty when practicing airway management. This versatile model ranges from challenging to very difficult to intubate, facilitating a variety of training scenarios from beginner to advanced users. The AirSim Difficult Airway has been designed to improve management of the unanticipated difficult intubation in accordance with the Difficult Airway Society guidelines.

Model Features

Difficult Airway

The internal airway has an enlarged tongue, swollen posterior cartilages and an elongated epiglottis to provide a challenging airway management user experience.

Laryngospasm

Adjustable from normal to complete airway obstruction.

Sudden closure of the vocal cords can be simulated.

Upon successful user action, laryngospasm can be mechanically reversed to provide positive feedback

Displaced larynx

Ability to displace the larynx into a more anterior position to make visualisation of the vocal cords more challenging. Adjustable to provide various levels of difficulty.

Tongue Edema

The tongue can be swollen to simulate various degrees of tongue angioedema.

Receding / protruding Mandible

Various degrees of malocclusion, one mandibular variant and two maxillary variants. Overbite and underbite features simulated

Trismus

Ability to grossly restrict mouth opening to make insertion of airway devices more challenging

Neck lock / Head rotation restriction

Ability to simulate cervical arthropathy through restricting the range of motion of cervical spine Neck lock / Head rotation restriction

Simulate needle and surgical cricothyroidotomy and percutaneous tracheostomy.

Easy identification of sternal notch, tracheal rings & cricoid and laryngeal cartilages Interchangeable larynx allows quick and easy repeated use Wraparound replaceable neck skin rotates for 10-15 incisions