

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz Bedsore and Wound care Model: BT-CSCB

Order code: 4119.BTCSCB



Information about product price on demand

Parameters

Quantitative unit ks

Sada obsahuje 4 fáze proleženin

- Stage 1: Unbroken skin / Redness of the skin
- Stage 2: Partial-thickness skin or tissue Loss / Formation of blisters
- Stage 3: Full-thickness skin loss / Damage to skin
- Stage 4: Full-thickness tissue loss / Exposed muscle, tendon and bone. Life-treatening infection



Wound is included in all buttocks manikin

- Dressing practice according to the stage of bedsores
- Can be used alone by not attaching to other simulators

Specification & Components

4 Models: Stage 1, 2, 3, 4 (1 model for each)

1 Wearable belt 1 Storage case 1 user manual

Size: $30 \times 28 \times 8$ cm (1 buttocks)

Weight: about 0,8 kg