



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in  
Hradec Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: info@helago-cz.cz  
Web: http://www.helago-cz.cz

## LF00855 - Male Catheterization Simulator

Order code: **4107.LF00855**



Cena bez DPH

1.071,00 Eur

Price with VAT

1.295,91 Eur

Parameters

Quantitative unit

ks

With the Life/form® Male Catheterization Simulator, catheterization can be practiced just as on a real patient. A lubricated catheter can be inserted in the urethral orifice, passed through the urethra, and into the bladder. When the bladder is successfully entered, artificial urine (water) will flow from the catheter. The student will feel the normal restrictions caused by the mucosal folds, bulbous urethra, and the internal urethral sphincter, just prior to entrance into the bladder. The experience teaches proper positioning and movement of the penis to allow the catheter to pass easily with a minimum of discomfort to the patient. Just as in a real patient, the cuff must be completely deflated before the catheter is removed. Improper use of a Foley can result in damage to the simulator, just as when using a Foley catheter with a real patient. Includes a hard carrying case, lubricant, catheter, and teaching guide. Three-year warranty. Carton size: 17" x 12" x 10" (43 x 30 x 25 cm). Sh. wt. 12 lbs. (5.44 kg).

### Available Supplies

- LF01127 - Foley Catheter
- LF01128 - Foley Catheter (Pkg. of 10)

- LF00985 - Lubricant Kit
- W09919 - REN Cleaner
- SB06172 - Fluid Administration Set