

Web: http://www.helago-cz.cz

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz

PP00031 - Traction Splint Trainer Order code: 4107.PP00031



Parameters Quantitative unit

ks

Traction splint application requires the ability to regulate how much force is applied to a patient's leg. There has never been a way to monitor that, until now. The Traction Splint Trainer can be set to correspond to an individual's body weight. This allows the instructor to determine if the student is correctly applying the rule of thumb that states "10% of the body weight traction results in adequate relief for, but no further damage, to the patient." Adjust for body weights ranging from 45 to 135 kg. The electrical, push-button panel makes simplified selection easy. This AC electrical trainer comes completely assembled and requires attachment to a base to secure it during training sessions. After applying traction, the deformity of the thigh recedes and splint can be applied, including any standard ankle hitch.