


HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

Swallowing Mechanism Model





Information about product price on demand

Parameters

Quantitative unit ks

Weight: About 12kg

Case size: W70 x D39 x H40cm

Accessories: Transparent cover, Pillow, Bolus model with wire, Baby powder, Storage case

Learn the mechanism of swallowing, and practice practical way of preventing miss-swallowing!

This Swallowing Mechanism Model is a medical teaching equipment that was developed with the purpose of helping helpers learn the mechanism of swallowing and the proper method of preventing swallowing disorders (miss-swallowing) when giving assistance in eating. It should help prevent pneumonia caused by swallowing disorders which is often overlooked in caring for the elderly.

The tongue is soft and able to pull out. Can have swallowing drills, study the tongue and muscles surrounding the mouth, and practice absorption.

- Neck angle (head angle) can be bent about 15° forward and about 45° backward.
- Can comprehend the relation between neck angle (head angle) and movement of windpipe / esophagus.

Practice:

- · Study the mechanism of swallowing
- Learn what causes miss-swallowing
- Study proper posture when eating and its relation with bed angle
- Study neck angle (head angle) and its relation with miss-swallowing
- Practice emergency care when food gets stuck
- Swallowing drills and practice in giving mouth care
- Learn how to insert tube through nose, and the intermittent tube eating method
- Study internal mouth absorbency