


Web: http://www.helago-cz.cz

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz

Sakamoto Tracheotomy Trainer Order code: 4109.M172

ks



Training of both conventional and percutaneous tracheotomy and puncture/incision into the cricothyroid ligament is available. For conventional tracheotomy, different types of incision can be performed including longitudinal, transverse, crucial, U-shaped, and inverted U-shaped incision. Simulating a patient in a supine position with the neck extended, the SAKAMOTO Tracheotomy Trainer allows the user to determine the proper incision site while identifying arterial locations, and to observe intratracheal conditions from the head.

Features

1. More than one time incision is possible

You can reuse the trachea by turning it upside down or inside out, and the skin by sliding it to either side or turning it upside down.

2. Identification of arterial locations

You can determine the proper incision site while identifying arterial locations.

3. Observation of intratracheal conditions

You can observe intratracheal conditions.

Practice

Conventional tracheotomy Percutaneous tracheotomy Puncture/incision of the cricothyroid ligament