

Sinus Lift Practice Model

Order code: **4112.345090**



Cena bez DPH

165,20 Eur

Price with VAT

199,89 Eur

Parameters

Types of replacements

Zubní implantáty

Quantitative unit

ks

- Jaw model for sinus lift practice
- The bone thickness of the posterior area is 3mm on the left side and 9mm on the right side for lateral window and osteotome method practice
- Schneiderian membrane is made of thin material
- Gingiva: gingiva with periosteum for use in incision/suture
- Bone Material: duplicates type III (1-layered construction)
- Tooth Model: removable anatomic tooth model
- Other: Schneiderian membrane

Training Method:

- Incision
- Lateral window technique
- Osteotome
- Incision, gingival ablation
- Suture
- Drilling
- Immediate implant placement

- Implant insertion
- Surgical guide plate
- Alginate impression
- Fabrication of superstructure