



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in Hradec Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: [info@helago-cz.cz](mailto:info@helago-cz.cz)  
Web: <http://www.helago-cz.cz>

**Small Uterus**  
Order code: **4113.UT001**



Information about product price on demand

#### Parameters

Obstetrics, gynecology - filter

Quantitative unit

Common gynecological and obstetric actions

ks

The Gynaecological Trainer is a highly realistic simulation of the female reproductive organs. Made from a very life-like material, it allows the medical practitioner to practice everyday procedures such as bi-manual palpation, pap smears, and the use of the speculum.

The Gynaecological Trainer has a unique clip-on system that enables the user to quickly and easily change over to a different uterus of a different size or pathology.

The clip is a simple and effective design. It's flexibility allows it to simulate the movement of the uterus on the cardinal ligaments when pushing on the cervix during bi-manual palpation. To change over a uterus, simply turn the model over, unscrew the clip from the back, swap over the uterus, and re-attach the clip.

The Gynaecological Trainer comes with seven different uteri. They are, from the top and left to right, hollow, small, medium, large, ovarian cyst, fibroids and retroverted.

The highly realistic nature of the anus on the the Gynaecological Trainer makes it possible, particularly for the student, to feel the close proximity of the bowel to the birth canal.

The Gynaecological Trainer is made of a very soft, highly elastic, and very durable material.

It is easy to clean, and following simple maintenance instructions ensures the models will last a long time.

They are also easy to repair if they are accidentally damaged.