

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz LF03674 - Child CRiSis™ Starter Torso

Order code: 4107.LF03674



Cena bez DPH 3.450,00 Eur
Price with VAT 4.174,50 Eur

Parameters

Construction Adult torso

Quantitative unit ks

- Airway Management
- Arrhythmia Recognition

- Defibrillation
- ECG Monitoring (4 sites)
- Manual Carotid Pulse
- Pacing

The Life/form® Child CRiSis™ Starter Torso can be completely upgraded to a dramatic, full PALS training manikin. Present your students with the challenges they may face in the real world. This airway management trainer is perfect for practicing skills on pediatric patients and represents the size of an 8-year-old patient. Practice intubation, ventilation, suction, CPR, and jaw thrust techniques. Manual carotid pulse. Realistic anatomy and landmarks include teeth, tongue, oral and nasal pharynx, larynx, epiglottis, arytenoids, false and true vocal cords, trachea, lungs, esophagus, and stomach. The trainer allows you to practice oral, digital, and nasal intubation, as well as E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and KING System insertion. Separate lungs for auscultation. Inflatable stomach bladder indicates esophageal insufflation. With its slightly anterior position, the Life/form® Child CRiSis™ Starter Torso is ideal for ACLS, nursing, paramedic, and EMT training at every level. This trainer with defibrillation chest skin lets you practice defibrillation using standard manual, automatic, or semiautomatic external monitor defibrillators. An internal load box absorbs the full strength of every shock to protect students and equipment. It is also possible to monitor the manikin, just like a real patient, and it is compatible with all standard brands and types of defibrillators, monitors, and patient simulators. Features four ECG sites and two defibrillation sites. This torso is designed to be the foundation of the Life/form® Child CRiSis™ Manikin. Arms and legs can be added later to provide IV access, intraosseous infusions, and blood pressure monitoring. If you do not have an ECG simulator and want to practice ECG recognition, you may want to consider adding the Life/form® Interactive ECG Simulator to your training (LF03670). If advanced airway management is preferred, consider starting with the Life/form® Advanced Child CPR/Airway Management Torso with Defibrillation Features (LF03763). Pump spray lubricant and instruction manual are included.