



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in Hradec Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: [info@helago-cz.cz](mailto:info@helago-cz.cz)  
Web: <http://www.helago-cz.cz>

**Elara 11 D Line**  
Order code: **2106.ELARA11230D**



Information about product price on demand

#### Parameters

Types of sterilizers	Parní
Volume [l]	28
Inner dimensions (diameter × d) [mm]	280 × 430
Quantitative unit	ks

The Elara 11 is designed for sterilization of wrapped, porous and hollow loads. The chamber volume is large, external dimensions are small, making it especially suitable for clinics requiring a high turnover. The control system is smart and intuitive, easy to use and has a large memory for consistent and reliable performance.

#### Efficient

- Short cycle and high volume pre and post vacuum pump for perfect class B sterilization results
- 316L grade stainless steel chamber and door with electro-polish finish
- A built-in steam generator provides steam and facilitates a standby steam mode
- The door closes smoothly and locks automatically

#### Features

- Large 28.5 Liter chamber within a small foot print. A perfect fit for standard 60 cm shelves or counters
- Provided with 5 stainless steel trays and can also hold cassettes
- Water filling from the front economizes on office space by eliminating the need to keep the top of the autoclave accessible

#### Safety features

- A cycle will not start unless the door is properly locked
- An automated double locking safety device prevents the door from opening while the chamber is pressurized
- An automatic safety shut-off prevents overheating of the chamber
- The autoclave comes with safety pressure-release valves for the chamber and generator