

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

Intramuscular Injection Training Buttocks Model Order code: 4119.BTCSIMT



- CSIM-T Simulators feature realistic shape and soft material with realistic tactile feeling.
- The buttocks model is composed of epidermis layer, muscular layer and skeletal frame.
- Ventrogluteal and dorsogluteal site injection practice.
- The skin is made of soft silicon that leaves no needle marks.
- Buttocks skin self-replacing.
- Both side of buttocks can be used for practice.
- Drug injection(Only Water) and Automated drain.
- User can check that the injection area is correct or not through the LED / warning sound.
- User can select AC 110V/220V 60/50Hz power or Battery power.