




HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec
Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

SMA AW201 - SmartMan ALS Airway CPR Pro

Order code: **4118.SMAAW201**



Information about product price on demand

Parameters

Construction

Quantitative unit

Adult torso

ks

A high fidelity airway manikin that provides life-like interaction for air movement related to pressures and movement in the chest. This brings more reality to your training and improves how chest compressions and ventilations are performed.

Both BLS and ALS training

Includes BLS and ALS training Real Time Feedback for BLS and CPR with Advanced Airway with feedback for 10:1 Asynchronous Continuous Compressions with a ventilation every 6 Seconds.

Delta Tidal Flow

Delta Tidal Flow Interaction System is found ONLY in SmartMan. This technology provides accurate measurement and excellent physical feel of the interaction of air movement. The real time digital display shows the volume, rate and timing air in and out of the lungs as it is related to movement of the chest.

Intubate

Practice Intubating and receive a score on how quickly and effectively you are able to place the device. Practice when chest compressions are ongoing and then move into 10:1 Asynchronous Continuous Compressions with a ventilation every 6 Seconds)

NOTE: Control PC or tablet is not a part of delivery.