



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court
in Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: <http://www.helago-cz.cz>

1022816 - Simulated Patient Monitor with Debriefing & CPR Feedback - REALiTi Pro

Order code: **4127.1022816**



Information about product price on demand

Parameters

Quantitative unit

ks

The simulated patient monitor REALiTi Pro is a highly advanced vital sign simulation system developed by iSimulate. REALiTi Pro is the advanced version of the REALiTi360 Ecosystem and allows medical educators to run multiple scenarios from basic to sophisticated. It is really a mobile complete sim center in a bag. You are not confined to a particular location and can conduct your training and evaluation in-situ in an ambulance, a helicopter, a hospital setting or in a sim center or skills lab. You are only limited by your imagination.

With REALiTi Pro, users get to choose 5 simulated proprietary patient monitor screens so you can provide training on the vital signs monitors you use in real life in hospital and in pre-hospital settings. In a medical or nursing undergraduate environment, you could even consider adding additional monitors to prepare your students for multiple monitor types they are likely to encounter during their careers.

REALiTi Pro includes all the hardware and software to run detailed debriefings and a complete CPR module that you can use with any patient manikin. REALiTi Pro is also the perfect tool for distance-based learning.

REALiTi Pro by iSimulate is a smart, integrated and modular simulation ecosystem for more reality and flexibility in your simulation training. It is a complete in-situ patient monitor simulator with live video streaming, CPR feedback and simulated patient records.

REALiTi Plus supports the training of the following courses, provided the corresponding Simulated Screen is used:

- ALS - Advanced Life Support
- ATLS - Advanced Trauma Life Support
- ACLS - Advanced Cardiac Life Support
- APLS - Advanced Paediatric Life Support
- EPALS - European Paediatric Advanced Life Support
- BASIC - Basic Assessment and Support in Intensive Care
- ILS - Immediate Life Support
- PALS - Paediatric Advanced Life Support
- ITLS - International Trauma Life Support
- PHTLS - Pre Hospital Trauma Life Support
- PILS - Paediatric Immediate Life Support
- PTLS - Paediatric Trauma Life Support
- NLS - Newborn Life Support
- TNCC - Trauma Nursing Core Course

REALiTi Pro mimics real defibrillators, monitors and ventilators:

Working with leading defibrillator and monitor manufacturers, iSimulate have created hyper-realistic interfaces which look and work like the real patient monitors.

- Meticulous attention to detail
- Looks, sounds and functions like a real patient monitor
- Standalone ventilator simulator can be used independently
- Ventilator works in portrait and landscape modes

Additional patient monitor screens are constantly being developed and added by the iSimulate team. Please have a look at the list below to see all patient monitors/ vital signs monitors and defibrillators to choose from:

- corpuls1 Patient Monitor Screen Simulation
- corpuls3 Patient Monitor Screen Simulation
- corpuls® AED Defibrillator Screen Simulation
- CARESCAPETM B40 Patient Monitor Screen Simulation
- LIFEPAK® 1000 Patient Monitor Screen Simulation
- LIFEPAK® 15 Patient Monitor Screen Simulation
- LIFEPAK® 20 Patient Monitor Screen Simulation
- Medtronic Capnostream™ 35 Patient Monitor Screen Simulation
- Philips IntelliVue MX800 Patient Monitor Screen Simulation
- Philips HeartStart MRx Emergency Care Patient Monitor Screen Simulation
- Philips HeartStart MRx for Hospital Patient Monitor Screen Simulation
- Welch Allyn Connex® VSM 6000 Patient Monitor Screen Simulation
- Zoll EMV+® Critical Care Ventilator Screen Simulation
- Zoll® Propaq® MD Patient Monitor Screen Simulation
- Zoll® R Series® Patient Monitor Screen Simulation

- Zoll® X Series® Patient Monitor Screen Simulation
- Schiller DEFIGARD Touch 7 Patient Monitor Screen Simulation
- Schiller PHYSIOGARD Touch 7 Patient Monitor Screen Simulation

Additionally, the following generic (unbranded) monitors are included:

- Generic Screen - AED Screen Simulation
- Generic Screen - Defibrillator Screen Simulation
- Generic Screen - Monitor Screen Simulation
- Generic Screen - ICU View 1 Screen Simulation
- Generic Screen - ICU View 2 Screen Simulation

REALiTi Pro comes with an integrated wireless video module for debriefing and distance-based learning:

- The Video module is integrated into REALiTi Pro by iSimulate
- Live stream video wirelessly
- Record and review video footage
- Clear head-up display of vital signs
- Video is transferred to Control at the end of a scenario and can be exported with vital signs parameters overlaid

REALiTi Pro includes a CPR module for real-time CPR feedback:

- Versatile CPR modules
- Can be worn on the wrist, placed into a puck or inside a manikin
- Real-time CPR feedback
 - Rate
 - Depth
 - Recoil
- CPR report for review
- Activates CPR features on proprietary screens such as Zoll and Corpuls

REALiTi Pro offers millions of waveforms, from a single touch:

Users can add ectopics, arrhythmias and artifacts to any ECG waveform.

- PVC (Premature Ventricular Contraction), PAC (Premature Atrial Contraction) and PJC (Premature Junctional Contraction) can be added to any ECG waveform
- Sinus arrhythmia linked to respiratory rate
- Artifacts and electrical interference
- Includes Critical Care Waveforms - CVP (Central Venous Pressure), ICP (Intracranial Pressure) and PAP (Pulmonary Artery Pressure)
- Waveform preview shows trends

REALiTi Pro has a built in radiology and 12 lead ECG library:

In collaboration with Life In The Fast Lane (<https://litfl.com/>) REALiTi offers a complete library of resources.

- Curated radiology images
- Dynamic 12 lead ECG library
- Advanced ECG lead control

REALiTi Pro fosters engaging assessment for students and distance based learning:

Easy to use checklist interface for simple student assessment. Add the REALiTi Engage App to allow observers to comment, vote and view media.

- The Checklist interface is completely configurable and easy to use
- Observers can use the REALiTi Engage App on their mobile devices

- Export a complete report of the simulation for student records

REALiTi Pro includes a built in Training Management System to collect all your simulation data in one central location:

- Event log, CPR compressions, vital signs history, scoring, waveforms, video and all other data saved on Control
- Video is transferred to Control at the end of a scenario and can be exported with vital signs parameters overlaid
- Browse past scenarios and students
- Export scenarios reports as a PDF to share
- Export all data as CSV to share

Simulated patient records become real with REALiTi Pro:

With REALiTi Chart you have a simulated patient record to view radiology, 12 lead ECGs and laboratory results in real time.

The REALiTi 360 Patient Monitor App - now available with new features in the V9 Release:

The iSimulate team has just announced the release of the V9 Build of the REALiTi360 patient monitor simulator. With the update, users now have the following new features available:

- Critical Care Waveforms - CVP, ICP and PAP
- Multiple ECGs are now 12 leads
- Classification of imported media on Control with options: General, Ultrasound, Laryngoscopy, Bronchoscopy, Radiography
- COVID 19 scenario

When media is chosen to be shown on the monitor, depending on the category chosen on Control, the media will now be displayed inside a simulated device (i.e. Laryngoscopy images now show inside a Karl Storz CMAC device). Whether or not static HR tones or saturation tones are played is now also configurable per device. Additional features and updates include:

- The waveform preview will now show trends
- The critical care waveforms now show on a separate graph in the PDF report
- The Schiller Defigard Auto-Screenshot has been added as a setting
- 12 leads have been added for Atrial Tachycardia and Atrial Flutters
- Significant improvements to 12 lead ECGs across the board
- 3rd Degree Heart Block 12 lead ECG
- A custom UI for displaying the critical care waveforms has been added to the GE CareScape
- The QRS volume function has been added to Corpuls3
- Scale options for arterial lines and CVP have been added to the GE CareScape
- Scenario details are now sent to the Chart as documents when the scenario starts
- The CPR game is available

Delivery content of REALiTi Pro by iSimulate:

- 5 branded patient monitor screens of your choice
- 1 "Patient Monitor Simulation Bag"
- 1 Simulation Accessory Set with:
 - 3 leadwires ECG cable snap
 - 7 leadwires ECG cable snap
 - 10 leads ECG truck cable
 - SpO2 sensor
 - Adult NIBP cuff
 - Defibrillator/ AED connector
 - AED Defibrillator pad
 - Disposable adult ECG electrodes
 - EtCO2 sensor gas drier tube
- CPR Module software and hardware (with sensors for manikins and standardized patients)
- Video debriefing software and hardware (incl. 1 Video Tablet)
- 1 WiFi Router
- 1 Patient Monitor Tablet with REALiTi360 App

- 1 Instructor Tablet with hard case and REALITI360 App