



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: <http://www.helago-cz.cz>

Loopster basic, rotator
Order code: **1305.04015000**



Cena bez DPH

972,00 Eur

Price with VAT

1.176,12 Eur

Parameters

Device type

Test tube shaker

Manufacturer

IKA

Maximum speed

20 rpm

Rotational Speed

fixed 20

Shaking motion

Rotational

Timer

NO

Maximum load	2 kg
Quantitative unit	ks

Rotator for smooth but effective mixing with a fixed speed. Ideal for biological samples, particularly blood samples. Suitable for separating and mixing processes of powdered and liquid samples in Eppendorf tubes and Greiner tubes up to 50 ml. The easy removable attachments allow simple and quick cleaning of the unit of spilled samples.

- Freely adjustable rotary angle
- Different attachments allow a variety of applications
- Usable at 4 to 50 °C in incubators
- Easy operation by foil keypad
- Suitable for continuous operation
- 20 rpm (fixed)