


HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

## DW-HL398HC Ultra Low Freezer For Laboratory And Medical Twin-compressor Freezer Order code: 2203.DWHL398HC



Twin-compressor, if one is damaged, the other one could keep the temperture at -70? stably. The digital temperature display can show various parameters simultaneously, such as temperature inside the cabinet, power voltage and environmental temperature, indicating operating status clearly The high-precision microcomputer temperature control system The well-developed alarm system has multiple alarm functions: high temperature/low temperature alarm, door ajar alarm, power failure alarm, sensor failure alarm, low battery alarm, alarm of poor heat dissipation of condenser, system failure alarm, etc. The keyboard lock function and password protection function can prevent any

adjustment of operating parameters without permission The rotating handle facilitates door opening The adjustable layered structure is convenient for storage of various things