


HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

7150 - Simulátor EGD Order code: 4003.7150



Information about product price on demand

Parameters

Quantitative unit

ks

This high-quality practice model for esophago-gastro-duodenoscopy allows insertion and navigation of an endoscope from esophagus to duodenum.

Special soft materials reproduce the gastrointestinal tract lifelike, allowing training of endoscopy as on real humans. The internal structures look lifelike in the endoscope and also behave like real tissue. The tongue, larynx, esophagus, stomach and duodenum are shown. The mouth, esophagus and duodenum can be swapped out and replaced with optional higher difficulty modules. The basic model comes in difficulty level 1, levels 2 and 3 can be ordered as accessories.

The model has polyps in 4 stages as well as an ulcer to practice navigation and diagnosis.

The following exercises can be performed:

- Navigation training from mouth to duodenum.
- Clipping (if this option is used, the stomach must be changed after some exercises)
- Biopsy (if this option is used, the stomach must be changed after some exercises)
- Sling training (if the sling is tightened completely, the stomach must be changed after some exercises)