

HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec
Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

Hydro H11V Multiple Water Bath

Order code: **2001.L003069**



Cena bez DPH

938,00 Eur

Price with VAT

1.134,98 Eur

Parameters

Effective height [mm]	100
Quantitative unit	ks

The special baths for gentle vaporization work from columns, Erlenmeyer flasks or beakers are suitable for unattended continuous operation in the laboratory thanks to water level controllers and a water shortage protection. LAUDA Hydro vaporization water baths are available in five models, all equipped with a removable perforated cover with a multi-part ring set made of heat-resistant plastic.



Working temperature min.
25 °C



Working temperature max.
100 °C

Technical Attributes

Working temperature min.	25 °C
Working temperature max.	100 °C
Ambient temperature min.	10 °C
Ambient temperature max.	40 °C
Temperature stability	3 ±K
Application	Inside the bath
Heater power	1.5 kW
Power consumption	1.5 kW
Bath volume max.	10.5 L
Bath depth	100 mm
Dimensions (WxDxH) in mm	450x300x155