

Hydro H19V Multiple Water Bath

Order code: **2001.L003070**



Cena bez DPH

1.072,00 Eur

Price with VAT

1.297,12 Eur

Parameters

Effective height [mm]

100

Quantitative unit

ks

The special baths for gentle vaporization work from columns, Erlenmeyer flasks or beakers are suitable for unattended continuous operation in the laboratory thanks to water level controllers and a water shortage protection. LAUDA Hydro vaporization water baths are available in five models, all equipped with a removable perforated cover with a multi-part ring set made of heat-resistant plastic.



Working temperature min.

25 °C



Working temperature max.

100 °C

Technical Attributes

Working temperature min.

Hydro H 19 V

25 °C

Working temperature max.

100 °C

Ambient temperature min.

10 °C

Ambient temperature max.

40 °C

Temperature stability

3 ±K

Application

Inside the bath

Heater power

1.5 kW

Power consumption

1.5 kW

Bath volume max.

18.4 L

Bath depth

100 mm

Dimensions (WxDxH) in mm

690x300x155