

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz LF00958 - Pediatric Arm Order code: 4107.LF00958



Cena bez DPH 714,00 Eur

Price with VAT 863,94 Eur

Parameters

Cannula, injection, puncturing - filter Horní končetiny

Quantitative unit ks

- · Practice Venipuncture and Intramuscular Injection Techniques and Procedures Used in Young Children
- Exact Reproduction of the Arm of a Six-Year-Old Child

A soft foam is used to simulate the deltoid muscle and allows the student to "get the feel" of administering intramuscular injections to children. Simulated bone in the shoulder defines and limits the injection area. Water maybe used as an injection fluid in the shoulder area. Includes a Life/form® Arm with white skin, fluid supply bag, 3 cc syringe, one pint Life/form® blood, teaching guide, and hard carrying case. Three-year warranty. Carton size: 22" x 11" x 6" (56 x 28 x 15 cm). Sh. wt. 10 lbs. (4.54 kg).

Available Supplies

- LF00986 Replacement Skin and Vein Kits
- LF01056 Vein Replacements Only
- LF00845 Venous Blood Quart. Sh. wt. 0.38 lbs. (0.17 kg)
- LF00846 Venous Blood Gallon. Sh. wt. 0.88 lbs. (0.40 kg)