

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz



Cena bez DPH 1.485,00 Eur

Price with VAT 1.796,85 Eur

Parameters

Description Pro 24 kyvet

Quantitative unit ks

In addition, you can program the CR 3200 to carry out 8 of your individual digestions at freely selectable temperatures up to 338 °F (170 °C).

**CR 3200** 

Order code: 2109.1P221

## **Optional accessories:**

• Quality assurance is constantly increasing in importance, even in the operational analysis sector. The CR 3200 and CR 4200 thermoreactors are both equipped with the external temperature sensor TFK CR (Order No. 250 100) as a testing aid. This temperature sensor can be plugged into the interface in place of a cuvette and the set and actual temperatures can be outputted either to a printer or a PC. This means that the function can not only be monitored, but also documented.