


HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

PCR box PW-01 with HEPA filter

Order code: 9002.AAHB3002



Cena bez DPH 3.135,00 Eur
Price with VAT 3.793,35 Eur

Parameters

Inner dimensions (w × d × h) [mm] $700 \times 585 \times 602$ [mm] Ext. dimensions (w×d×h) $720 \times 605 \times 964$ [mm]

Quantitative unit pcs

Using a HEPA filter provides an effective particulate-free working space for ambient environment conditions.

Effective sterilization 254nm ultraviolet lamp and Fluorescent lamp

- Adjustable UV-timer for precise control of decontamination
 - o The variable digital timer for UV exposure with detail settings by second. (Delay OFF : Max. 30min)
 - $\circ~$ Long life (Up to 8,000 hour) UV lamp
- Anti-glare fluorescent lamp is installed on the center ceiling of the work space area, minimizing shadows and relieves eye strain

User-Friendly Control System

- Microprocessor-based control system
- Bright and vivid LCD display with responsive touch buttons
- Easily indicating the total use time and operating time on the screen

A reliable interlock safety system

- An audible alarm for product and sample protection against UV over exposure
- For safety UV-opaque and resistant 99,99%, thick acryl front, back and sides
- Automatic UV lamp safety shut OFF, in case the door is opened, which prevents accidental UV exposure
- Self-protection against over-current