



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in Hradec Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: [info@helago-cz.cz](mailto:info@helago-cz.cz)  
Web: <http://www.helago-cz.cz>

## Defibrillator iPAD CU-SP2

Order code: **4401.P04986**



Information about product price on demand

### Parameters

Types of defibrillators and AED

Ostré defibrilátory

Quantitative unit

ks

Defibrillator iPAD CU-SP2 is easy to use and clear defibrillator with colored display. It's maximally effective and life saving. In semi-automatic mode the defibrillator guides user step by step with simple and clear orders in Czech language. In manual mode the user assesses whether the patient needs a defibrillation shock or not. CU-SP2 allows for a synchronization of shock with R wave. It's light and powered by battery, which allows for the maximal mobility.

Defibrillator i-PAD CU-SP2 is intended for use during ventricular fibrillation and quick ventricular tachycardia. Thanks to an automatic impulse volume adaptation the defibrillator is suitable for use in noisy environment.

You can order a carry case for this defibrillator (4401.P04415).