



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in Hradec Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: [info@helago-cz.cz](mailto:info@helago-cz.cz)  
Web: <http://www.helago-cz.cz>

**ANGIO Mentor**  
Order code: **4506.ANGIO**



Information about product price on demand

## ANGIO MENTOR PLATFORMS

The ANGIO Mentor family of products is an excellent example of our commitment to provide flexible, cost-effective solutions for a wide range of settings.

Any software version we release is compatible with all our platforms, including the Procedure Rehearsal Studio.

### ANGIO Mentor Suite

Provides a true-to-life Cath Lab/ ANGIO Suite clinical environment with a full body mannequin and a large fluoroscopy/vital signs monitor – ideal for multi-disciplinary team training.

The system allows tool introduction from 5 optional access sites, and is ergonomically designed to allow table height adjustment and C-arm control panel position.

Integration of an actual C-arm into the simulation environment is an optional addition to help create an even more realistic setting and to enhance team training for endovascular procedures. In this set-up, the physical rotation of the C-arm is reflected in both the fluoroscopic image and the virtual C-arm simulation on the monitor.

## **TEE Probe Integration- Coming Soon**

Creates a realistic work environment and ideal for team training!

A modified TEE probe is connected to adapted patient mannequin, which allows the trainee to introduce the TEE probe and practice basic TEE operation, maneuvering and image interpretation. Provides practice of live echo guidance for interventions that require TEE imaging.

### **TEE probe kit includes:**

- Modified mannequin for TEE practice
- Probe and related hardware
- Specific SW package:
  - TEE training module
  - Revised Transseptal Puncture module with live TEE functionality

## **ANGIO Mentor Flex**

A compact and easy to use solution for various training needs, especially suitable for remote meetings and courses. Light, foldable and packed in a small case, it can be checked-in as luggage on flights within airline weight restrictions.

The system is easily set up by connecting only 3 cables. Smooth service and field support are enabled by the dockable tracking stations that can be replaced quickly and easily without opening the system.

Available in both single and dual access configurations, the ANGIO Mentor Flex can run all ANGIO Mentor modules and the PROcedure Rehearsal Studio™ software.

## **ANGIO Tab Pro**

A light weight, carry-on solution, packed in a small and light-weight suitcase designed for maximum protection.

Ideal for device demonstration (for clinical reps), teaching procedural steps (for hospitals/training centers) or planning and rehearsing an upcoming procedure using a dedicated simulation file created with Procedure Rehearsal Studio software.

### **Runs complete procedure simulation:**

Enables physical introduction and manipulation of one tool using a tracking unit, as well as device deployment using an authentic handle. Other interventional devices are virtually introduced simultaneously controlled by keyboard and touch gesture

