

BT-CSCF - Female Catheterization Simulator

Order code: **4119.BTCSCF**



Cena bez DPH

1.936,00 Eur

Price with VAT

2.342,56 Eur

Parameters

Quantitative unit

ks

Training environment more similar to clinical situation

- We make a training environment which is more similar to clinical situation by using a wide range of catheters and a water-soluble lubricant jelly
- Realistic resistance during insertion of urinary catheter
- When a catheter enters the bladder, artificial urine will come out

- External genitalia is separable from the buttocks/Easy to maintain
- Indwelling Catheterization and Intermittent Catheterization
- Practice with various size(Fr) of catheter
- Dorsal recumbent position; the General position for a catheterization
- Warning sound when a catheter inserted into vagina
- Detachable extension legs

Features

- Life-like practical training
- Make position similar to the real patient
- Any water-soluble Jelly being used at clinical situation can be applied
- Training with various size of catheter is available
- Realistic practice with attachable extended legs
- Non-return valve : Prevent leaking after taking catheter out.
- Checking fluid after inserting a catheter into bladder correctly
- Alarm when a catheter inserted into a vagina ; to operate, a battery is needed
- Genital organs and buttocks are attachable : easy to maintain

Specification

- Adult buttocks model: 3.5kg (430×280×200mm)
- Extended legs (Left/Right): Each 1kg (320×150×150mm)
- Storage case: 2.4kg (440×330×390mm)
- Total: Approximately 8kg