


HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

M93C - CVC Insertion Simulator

Order code: 4103.M93C



Cena bez DPH

Price with VAT

2.268,00 Eur

2.744,28 Eur

Parameters

Cannula, injection, puncturing - filter

Quantitative unit

Catheter insertion

ks

Features

- Insertion practice is repeatable to obtain the best technique for patient safety
- Both landmark guided insertion and ultrasound guided techniques can be practiced
- The frictionless tissue-like simulation pad to facilitate the realistic sensation of CVC placement

Skills

- Perform hand hygiene practices and prepare aseptic field
- Decide on insertion site via ultrasound guidance or observation and palpation of landmarks
- Insertion under ultrasound guidance or via landmark technique
- Confirm flashback
- Advance guidewire and sheath
- Insert catheter

Set includes

- 1 male upper torso manikin
- 2 CVC placement pad
- 1 skin for cannulation training
- 1 transparent anatomical block
- 2 artery tube (transparent): 2 type
- 2 vein tube (blue): 2 type
- 2 vein pipe
- 1 pulsation unit (for landmark technique)
- 1 irrigation bottle (for UG technique)
- 1 plastic jar
- 1 sample needle
- 1 syringe (2.5ml, 50ml each)
- 1 introductory ultrasound training block