

61104 - Episiotomy Trainer Module

Order code: **4104.61104**



Cena bez DPH

337,00 Eur

Price with VAT

407,77 Eur

Parameters

Quantitative unit

ks

OVERVIEW

- Baby head has realistic look and feel, providing landmarks and appropriate tactile response when inserting fingers to guard the baby's head
- Repair tape increases the longevity of the episiotomy pad, allowing for multiple uses per pad
- Illustration pack to help trainees get familiarized with key landmarks and where to place the cut

REALISM

- Realistic representation of tissue to represent a stretched perineum

CLEANING

- Skin surface is washable using soap & water

SAFETY

- Product is latex free

SKILLS GAINED

- Identification of fontanelles
- Handling of a stretched perineum
- Reducing the risk of damaging the baby or mother
- Infiltration of perineum prior to episiotomy
- Performance of mediolateral or midline episiotomy
- Use of instruments

The Perineal Repair Stand (61109) in which the module fits, must be purchased separately.