

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Fax: 495 220 154 GSM gate: 602 123 096

Phone: 495 220 229, 495 220 394

E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

S415.250 Blood Pressure Reading Skills Trainer with OMNI® 2

Order code: 4108.S415250



Cena bez DPH 2.731,00 Eur Price with VAT 3.304,51 Eur

Parameters

Quantitative unit ks

The Blood Pressure Training System with OMNI® 2 includes a full-size adult left arm that may also be attached to select Gaumard® adult simulators. This is a versatile training tool developed to assist health professionals teach the processes and skills required to perform blood pressure auscultation procedures and techniques. The package includes the new OMNI® 2 control tablet for simple, tap-and-go operation.

The Blood Pressure Training System with OMNI® 2 is designed for training the processes and skills required to perform blood pressure auscultation procedures and techniques. Now control the arm wirelessly with just a tap using the OMNI® 2.

Automatic Korotkoff Sounds

Auscultate Phase 1 - Phase 5 Korotkoff sounds when cuff pressure is between programmed systolic and diastolic pressures with real stethoscopes

Automatic Radial Pulse

Programmable, radial pulse synchronized with heart rate is palpable when cuff pressure is less than the selected systolic blood pressure

Auscultatory Gap

Programmable auscultation gap automatically silences Korotkoff sounds

Speaker Option

Speakers with volume control allow students to hear the same Korotkoff sounds the individual student hears while auscultating with the stethoscope

Training with OMNI® 2

The new OMNI® 2 features simulation tools designed to help boost skill acquisition, enhance the simulated experience, and support assessment and debriefing. It's the one handy tool that can do it all.

OMNI® 2 Wireless Controller

Connect to the Blood Pressure Training System from up to 12 meters using the Bluetooth connectivity feature of the OMNI® 2 Handheld Tablet

Intuitively Easy Control

Change the systolic pressure, diastolic pressure, heart rate, and auscultation gap with just a tap