


## **AC6/6 - Bladder Phantom for Endoscopic Examination**

Order code: **4101.AC6/6**



Information about product price on demand

### Parameters

Quantitative unit

ks

Phantom comprises:

- Natural size trunk made of special plastic
- On a base
- External female genitalia with urethra and stretchable opening to the bladder (removable inspection window and closure cap)
- Six interchangeable parts demonstrating characteristic bladder diseases:
  - Papillary bladder tumour
  - Trabeculated bladder
  - Urinary calculus
  - Broad-base bladder tumour
  - Fibrinous cystitis
  - Radiation cystitis
- For examination, the endoscope is inserted through the urethra via a stretchable opening in the bladder
- By changing the lens – pro-grade urethra optic or bladder optic – the shaft can remain in the urethra and bladder
- Examination routine can be checked through the inspection window

### Optional Accessories

- AC 6/8 - Pathological Changes